

Buffer Solution
Baker Analyzed
pH 7.00 [Phosphate] Green



Material No: 5647.0500
Effective from: 2016/02/15

Product Specification

Test	Specification
pH at 20°C	6.98 – 7.02

For Laboratory, Research or Manufacturing use.
This buffer is traceable to NIST Standard Reference Buffer. Contains Potassium Dihydrogen Phosphate / Disodium Hydrogen Phosphate. Contains 0.01% Sodium Benzoate.

A handwritten signature in blue ink, appearing to read 'B. Rykowska'.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000