

Contec® Products for Cleaning Small Environments

Contec® Cleanroom provides a full solution product portfolio for cleaning and disinfecting biosafety cabinets, RABS, laminar flow hoods and isolators. When utilized in conjunction, these products:

- Prevent cross-contamination from operator over-reaching into the environment
- Promote best practice with overlapping strokes for even coverage
- Enforce compliance with Annex 1 to promote cleaning, disinfection and residue removal
- Enhance ergonomics due to lightweight tool and handle



EasyReach™ Cleaning Tool

- Lightweight for single handed operation and reducing operator fatigue
- Rectangular head allows for overlapping strokes
- Dry and presaturated polyester pads with good abrasion resistance and chemical compatibility
- Versatile to use with various QuickConnect™ handles to reduce strain and over-reaching

Part No.	Description		Size	Packaging
Hardware				
76417-122	EasyReach Stainless Steel One-Piece Lay Flat Frame		16.5" (42 cm)	1/case
76417-120	EasyReach Stainless Steel Lay Flat Mop Frame		7.5" (19 cm)	1/case
76480-378	QuickConnect EasyReach Polypropylene Lay Flat Frame		7.5" (19 cm)	1/case
77786-120	QuickConnect Carbon Fiber Handle		18" (46 cm)	1/case
77322-292	QuickConnect Carbon Fiber Telescopic Handle		18" - 32" (46-81 cm)	1/case
Cleaning Pads				
76417-134	EasyReach Sterile Ultrasonically Sealed Polyester Cleaning Pads	 Sterile	7.5" x 2.8" (19 x 7 cm)	4/bag; 25 bags/case
76449-778	EasyReach Sterile Cleaning Pads Quilted polyester	 Sterile	7.5" x 2.8" (19 x 7 cm)	4/bag; 15 bags/case
89094-582	EasyReach Sterile 70% IPA Presaturated Polyester Cleaning Pads	 Sterile	7.75" x 2.75" (20 x 7 cm)	8/bag; 8 bags/case

Disinfectants & Solutions

- Daily and sporicidal products for enhanced contamination control
- Cleaning with IPA reduces and removes residue build-up from surfaces



Part No.	Description		Size	Packaging
Intermediate Disinfectant				
76382-888	Contec Sterile CyQuanol Quaternary Ammonium/Ethyl Alcohol Cleaner Disinfectant		32 oz (0.9 L)	12/case
Sporicidal Disinfectant				
10148-324	Contec Sterile PeridoxRTU® Blend of hydrogen peroxide and peroxyacetic acid		32 oz (0.9 L)	6/case
Clean and/or Rinse				
89097-626	Contec® Sterile 70% IPA		16 oz (0.47 L)	12/case

Wipes for Irregular Surfaces

- Wipe down valves, bars and other irregular surfaces
- Dry wipes can be used with IPA or disinfectants
- Presaturated IPA wipes offer controlled release of solvent to surfaces



Part No.	Description		Size	Packaging
33503-128	PROSAT® Sterile™ Polynit Heatseal Wipes Presaturated with 70% IPA and 30% DI Water		9" x 9" (23 x 23 cm)	20/pouch; 40 pouches/case
68310-183	Sterile Polynit Heatseal Wipes Flat Stacked		9" x 9" (23 x 23 cm)	25/bag; 48 bags/case


avantorsciences.com
Any prices, product, and/or services details are current when published and subject to change. Offers subject to product availability. Void where prohibited by law or company policy. Any claims concerning products manufactured by third parties are the sole claims of the manufacturer and not those of Avantor or its affiliates. All products are subject to the terms and conditions of sale of Avantor or its affiliates, including the applicable limitations of liability contained therein. The information herein should not be considered a warranty of any kind. Please carefully review any product information sheets and safety information for the product. Trademarks are owned by Avantor Inc. or its affiliates, unless otherwise noted. | Visit [avantorsciences.com](https://www.avantorsciences.com) to view our privacy policy and additional disclaimers. © 2025 Avantor Inc. All rights reserved.