

Certificate of Conformity

VWR Biohazard Bag, 12 L Red with Printed Sign and Temperature Indicator

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No: 129-0033

Lot No: 0000954070906

Material: Polypropylene

Characteristics: Autoclavable upto 135°C

Quality System Compliance:

This certifies that the above listed product has been inspected and released in accordance with Plant Standard Production Methods and Process Quality Assurance Specifications.

QC Testing:

Established physical properties are tested in accordance with internationally recognized ASTM Test methods for Dart Impact Strength (ASTM D1709) and tear resistance (ASTM D1922).

Product Specifications:

The product allows use for autoclave sterilizer operating temperatures of 274 °F (134 °C).

BSE/TSE Statement: The film and resin suppliers confirm that some films, resins and masterbatch additive packages used to manufacture SUPROP® type film products and other films, may contain byproducts of bovine sourced raw materials in their manufacture. As advised by these suppliers, these animal byproduct additives are processed under highly refined and controlled industrial chemical processes and at very high temperatures prior to purification and subsequent use as a polymer additive.

This applies to all products comprised of either polypropylene or low-density polyethylene polymers. The animal derived byproducts are not added during the process of manufacture of the product but originate either in the manufacture of the resins or in the additives used in film formulations.

Certificate of Conformity

Latex Free Statement: The product is Latex Free

REACH Compliant: Yes

Manufacturing Date (DD/MM/YYYY): 11/04/2020

This document has been produced electronically and is valid without a signature

Date of issue: **April 2021**