

PRODUCT SPECIFICATION/CERTIFICATE OF CONFORMITY

PRODUCT CODE	3295/12AS
DESCRIPTION	10ML BURETTE CLASS AS WITH PTFE KEY
MANUFACTURER	PYREX®
MATERIAL	BOROSILICATE GLASS TYPE 3.3 COMPLIES WITH ISO3585 EUROPEAN AND US PHARMACOPOEIA GLASS TYPE 1
COMPLIANCE	PRODUCT COMPLIES WITH ISO385
PRINTED	SCREEN PRINTED IN WHITE ENAMEL
GRADUATIONS	0.02ML
TOLERANCE	± 0.02ML
OVERALL LENGTH	820MM
PACK QUANTITY1	1

This is to certify that the goods detailed above have been inspected and tested in accordance with our BS EN ISO9001 and BS EN ISO13485 quality systems, and conform in all respects to the acceptable quality levels (AQL) detailed therein and any relevant specifications or drawings, unless recorded below as concessionary goods.

Signed: P Burton
Quality Technician

Date: 24 June 2021