

# Certificate of Conformity

|                                         |                                                                                                             |
|-----------------------------------------|-------------------------------------------------------------------------------------------------------------|
| <b>VWR Catalogue Number</b>             | 113-1463                                                                                                    |
| <b>Description</b>                      | VWR® Basic Beard covers, White, Elastic Length 15.2cm, Stretched Length 40.6cm and width 25.4cm, with Latex |
| <b>Lot Number</b>                       | A5221747                                                                                                    |
| <b>Date of Manufacture (yyyy-mm-dd)</b> | 2022-11-07                                                                                                  |
| <b>Country of Origin</b>                | Manufactured in China                                                                                       |
| <b>Date of Issue (yyyy-mm-dd)</b>       | 2023-08-01                                                                                                  |

## Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

|                              | <b>Test Method</b> | <b>MD</b> | <b>CD</b> |
|------------------------------|--------------------|-----------|-----------|
| Break Tensile Properties [N] | ASTM D5034         | 27.61     | 17.97     |
| Elongation rate [%]          | ASTM D5034         | 73.28     | 68.81     |

## Product Specifications

|                               |                                                             |
|-------------------------------|-------------------------------------------------------------|
| <b>Material</b>               | Spunbonded Polypropylene Fabric (SPP) 10 g/m <sup>2</sup> . |
| <b>Sterilization</b>          | Not Applicable.                                             |
| <b>ATP Assay</b>              | Not Applicable.                                             |
| <b>DNase &amp; RNase Free</b> | Not Applicable.                                             |

|                                |                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE). |
| <b>Cytotoxicity</b>            | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Non-Pyrogenic Statement</b> | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Latex Statement</b>         | Product contains Latex.                                                                                                                                                                                                                                                                                                                                     |
| <b>BPA Statement</b>           | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>DEHP Statement</b>          | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>RoHS</b>                    | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>REACH Statement</b>         | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage Conditions</b>      | Room temperature, dry and shady.                                                                                                                                                                                                                                                                                                                            |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.