

Calbiochem[®]

Certificate of Analysis

Transferrin, Holo, Human Plasma, Low Endotoxin

Batch Number:	D00127502
Material Number:	616424-100MG
Molecular Weight:	80000
CAS Number:	11096-37-0

Analytical Data

Test	Tolerance	Result
Formulation:		Lyophilized in 20 mM sodium phosphate, pH 7.4, 150 mM NaCl
Solubility:		H ₂ O. Please refer to the vial label for the pre-lyophilization concentration.
Form:		Solid
SDS PAGE purity:	≥95.0 %	>98.0 %
"HBsAg; HIV,HCV antibodies":		Negative
Concentration before lyophilization:		119.20 mg/mL
Endotoxins:	≤0.1 EnU/mg	<0.1 EnU/mg

Storage and Handling: -20°C



Tofazzal Hossain, Quality Assurance Manager

12-Jun-2012
Date

USA and Canada
Tel (800) 628-8470
technical@calbiochem.com

France
Freephone
0800 126 461

Germany
Freecall
0800 100 3496

Europe
Ireland
Toll Free
1800 409 445

United Kingdom
Freephone
0800 622 935

All Other
European Countries
+44 115 943 0840

All Other Countries
Contact Your Local
Distributor
www.calbiochem.com
technical@calbiochem.com

techservice@merckbio.eu

A Brand of EMD Chemicals Inc., an affiliate of Merck KGaA, Darmstadt, Germany
www.calbiochem.com

FOR RESEARCH USE ONLY. NOT FOR HUMAN OR DIAGNOSTIC USE.