

Instrument basket with drop handles

P/N: 7500053 | SiSch 240/250/50 SR FG DM

HUPFER
we make work flow

Technical data



Frame wire thickness:	5.0 mm
Drop handle wire thickness:	4.0 mm
Side and base mesh size:	4.64 × 4.64 × 1.0 mm
Weight:	0.54 kg
Width:	250 mm
Depth:	240 mm
Height:	50 mm

Similar to illustration, technical modifications reserved. Without decoration.

Strainer bowl for transport, storage and provision of various sterile supply instruments.

Bowl in hygienic design for easy cleaning, made of high-quality electro-polished stainless steel round material. Special shape of base frame allows stacking of bowls onto one another. Easy handling thanks to two interior chest handles at the face sides.

Time and date of the request:
13.11.2020, 10:03:00

All information / dimensions are approximate, technical changes reserved. © Hupfer