

Certificate of Compliance

Catalogue Number: 189171

Lot Number:

1	6	9	0	4	0	8
---	---	---	---	---	---	---

Date of Manufacture 02 April 2021

Use by Date:- 01/04/2026

Product description: 7ml Container-Ps-No Label, Plastic Cap
STERILE

Component Materials: -

Container -Polystyrene
Cap -Polypropylene

Sterilisation: -

The lot has been produced under aseptic conditions. The aseptic conditions are achieved by manufacture into an ISO class 7 clean room conforming to BS EN ISO 14644-1 (1999) The lot has been sterilized by a post-production ethylene oxide process

This manufacturing lot has been sampled, physically examined tested, packed, stored and transported in accordance with Standard Operating Procedures as contained in the Greiner Bio-One quality manual relating to our ISO9001:2015 registration.



Noel Ennett
Operations Director