

PDS No. 1122xx	PRODUCT DATA SHEET	Page 1 of 1
Revision 04	Tube, 4 ml, PP, Round Bottom	 greiner bio-one
	Item-No. 112201	

1.	Description / Specification	
1.1	Description	PP tube with round bottom, 4 ml Appropriate grip stopper, Cat.-No. 303321
1.2	Dimensions	See customer drawing
1.3	Volume	<u>Working volume:</u> 3 ml
1.4	Material / Resin	PP (Polypropylene), free of heavy metal
1.5	Colour	translucent
1.6	Sterilization	No
1.7	Quality Control	<u>Product-Control:</u> testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-196°C to +121°C
2.3	Autoclavability	121°C / 20 min.
2.4	Centrifugation, max. RCF	5800 x g: swinging-bucket rotor 30000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information→Chemical Resistance of Resins
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	240
3.2	Pieces / Box	3600
3.3	Lot-No.	A YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information	
	-	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 04	Date 31 July 2014	Date 1 August 2014	Date 1 August 2014	
Date 18.10.2013	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Schulz	