

32 gauge, Small Hub RN Needle, 2 in, point style 3, 6/PK

Specification Sheet (Part/REF # 7762-05)



Small Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 100 µL and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

Product Specifications

Gauge	32 gauge
Point Style	3
Needle Length	2.00 in (50.8 mm)
Autoclavable	Yes
Gas Sterilizable	Yes
Sales Quantity	6 needles per pack
Hub	Small RN

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/needles/7762-05>

Spec. Version B