

Axygen® Universal Pipette Tips

Redesigned for the modern researcher.

For more than 20 years, Axygen Universal pipette tips have set the standard for quality, accuracy, and precision in laboratories worldwide. Now redesigned to meet every researchers needs, the next generation of Axygen Universal pipette tips deliver improved pipettor compatibility, durable hinged racks, and a complete offering from 10 µL to 1200 µL. Featuring tips with extended length, wide bore designs for viscous solutions, and our proprietary Axygen MaxRecovery+™ technology, Axygen Universal pipette tips are sure to meet and exceed your liquid handling expectations.

- ▶ RNase-/DNase-free
- ▶ Human gDNA-free
- ▶ Nonpyrogenic
- ▶ USP Class VI polypropylene

Axygen Universal Filter Pipette Tips

Racked, Filtered

	Standard Surface VWR Corning Cat. No.	MaxRecovery+ VWR Corning Cat. No.	Volume (µL)	Pack Type	Sterile	Quantity
	78692-602 PTF-300-R-S	78692-600 PTF-300-L-R-S	0.5 - 10	Fine point	Yes	96/rack, 10 racks/pack, 5 packs/case
	78692-622 PTXLF-10-R-S	78692-620 PTXLF-10-L-R-S	0.5 - 10	Extended length, graduated, fine point	Yes	96/rack, 10 racks/pack, 5 packs/case
	78692-598 PTF-20-R-S	78692-596 PTF-20-L-R-S	1 - 20	Beveled	Yes	96/rack, 10 racks/pack, 5 packs/case
	78692-584 PTF-100-R-S	78692-582 PTF-100-L-R-S	1 - 100	Beveled	Yes	96/rack, 10 racks/pack, 5 packs/case
	78692-592 PTF-200-R-S	78692-590 PTF-200-L-R-S	1 - 200	Graduated, fine point	Yes	96/rack, 10 racks/pack, 5 packs/case
	--	78692-594 PTF-205-WB-L-R-S	1 - 200	Wide bore	Yes	96/rack, 10 racks/pack, 5 packs/case
	78692-606 PTF-350-R-S	78692-604 PTF-350-L-R-S	1 - 300	Graduated, fine point	Yes	96/rack, 10 racks/pack, 5 packs/case
	78692-578 PTF-1000-R-S	78692-576 PTF-1000-L-R-S	100 - 1000	Beveled, graduated	Yes	96/rack, 10 racks/pack, 5 packs/case
	--	78692-580 PTF-1005-WB-L-R-S	100 - 1000	Wide bore	Yes	96/rack, 10 racks/pack, 5 packs/case
	78692-586 PTF-1200-C-R-S	78692-588 PTF-1200-L-R-S	100 - 1200	Fine point	Yes	96/rack, 8 racks/pack, 4 packs/case

Reload, Filtered

	78692-452 PRFLF-300-C-S	78692-450 PRFLF-300-C-L-S	0.5 - 10	Fine point	Yes	96/deck, 5 decks/pack, 10 packs/case
	78692-466 PRFLTXLf-10-C-S	78692-464 PRFLTXLf-10-C-L-S	0.5 - 10	Extended length, graduated, fine point	Yes	96/deck, 5 decks/pack, 10 packs/case
	78692-448 PRFLF-200-C-S	78692-446 PRFLF-200-C-L-S	1 - 200	Graduated, fine point	Yes	96/deck, 5 decks/pack, 10 packs/case
	78692-444 PRFLF-1000-C-S	78692-442 PRFLF-1000-C-L-S	100 - 1000	Beveled, graduated	Yes	96/deck, 5 decks/pack, 10 packs/case

Axygen® Universal Pipette Tips

10 µL Fine Point



Standard Surface VWR Corning Cat. No.	MaxRecovery+™ VWR Corning Cat. No.	Volume (µL)	Pack Type	Sterile	Quantity
78692-552 PT-300	78692-554 PT-300-L	0.5 - 10	Bulk	No	1,000/pack, 10 packs/case
78692-560 PT-300-R	78692-556 PT-300-L-R	0.5 - 10	Racked	No	96/rack, 10 racks/pack, 5 packs/case
78692-562 PT-300-R-S	78692-558 PT-300-L-R-S	0.5 - 10	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
78692-438 PRFL-300-C	78692-440 PRFL-300-C-L	0.5 - 10	Reload	No	96 tips/deck, 10 decks/pack, 5 packs/case

10 µL Extended Length, Graduated, Fine Point



78692-608 PTXL-10	78692-610 PTXL-10-L	0.5 - 20	Bulk	No	1,000/pack, 10 packs/case
78692-616 PTXL-10-R	78692-612 PTXL-10-L-R	0.5 - 20	Racked	No	96/rack, 10 racks/pack, 5 packs/case
78692-618 PTXL-10-R-S	78692-614 PTXL-10-L-R-S	0.5 - 20	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
78692-460 PRFL-TXL-10-C	78692-462 PRFL-TXL-10-C-L	0.5 - 20	Reload	No	96/rack, 10 racks/pack, 5 packs/case

200 µL Clear



78692-510 PT-200-C	78692-512 PT-200-C-L	1 - 200	Bulk	No	1,000/pack, 10 packs/case
78692-522 PT-200-C-R	78692-514 PT-200-C-L-R	1 - 200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
78692-524 PT-200-C-R-S	78692-516 PT-200-C-L-R-S	1 - 200	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
78692-526 PT-200-C-STK	78692-518 PT-200-C-L-STK	1 - 200	Stack rack	No	200/deck, 5 decks/ pack, 10 packs/case
78692-528 PT-200-C-STK-S	78692-520 PT-200-C-L-STK-S	1 - 200	Stack rack	Yes	200/deck, 5 decks/ pack, 10 packs/case
78692-434 PRFL-200-C	78692-436 PRFL-200-C-L	1 - 200	Reload	No	96 deck, 10 decks/pack, 5 packs/case

200 µL Beveled, Yellow



78692-530 PT-200-Y	78692-532 PT-200-Y-L	1 - 200	Bulk	No	1,000/pack, 10 packs/case
78692-538 PT-200-Y-R	78692-534 PT-200-Y-L-R	1 - 200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
78692-540 PT-200-Y-R-S	78692-536 PT-200-Y-L-R-S	1 - 200	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
78692-542 PT-200-Y-STK	78692-542 PT-200-Y-STK	1 - 200	Stack rack	No	200/deck, 5 decks/ pack, 10 packs/case

200 µL Wide Bore



--	78692-546 PT-205-WB-C-L	1-200	Bulk	No	1,000/pack, 10 packs/case
--	78692-548 PT-205-WB-L-R	1-200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
--	78692-550 PT-205-WB-L-R-S	1-200	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case

Axygen® Universal Pipette Tips

300 µL Graduated, Fine Point



Standard Surface VWR Corning Cat. No.	MaxRecovery+™ VWR Corning Cat. No.	Volume (µL)	Pack Type	Sterile	Quantity
78692-564 PT-350-C	78692-566 PT-350-C-L	1 - 300	Bulk	No	1,000/pack, 10 packs/case
78692-572 PT-350-C-R	78692-568 PT-350-C-L-R	1 - 300	Racked	No	96/rack, 10 racks/pack, 5 packs/case
78692-574 PT-350-C-R-S	78692-570 PT-350-C-L-R-S	1 - 300	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case

1000 µL Graduated, Beveled, Clear



78692-480 PT-1000-C	78692-482 PT-1000-C-L	100 - 1000	Bulk	No	1,000/pack, 10 packs/case
78692-488 PT-1000-C-R	78692-484 PT-1000-C-L-R	100 - 1000	Racked	No	96/rack, 10 racks/pack, 5 packs/case
78692-490 PT-1000-C-R-S	78692-486 PT-1000-C-L-R-S	100 - 1000	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
78692-430 PRFL-1000-C	78692-432 PRFL-1000-C-L	100 - 1000	Reload	No	96/deck, 10 decks/ pack, 5 packs/case

1000 µL Graduated, Beveled, Blue



78692-468 PT-1000-B	78692-470 PT-1000-B-L	100 - 1000	Bulk	No	1,000/pack, 10 packs/case
78692-476 PT-1000-B-R	78692-472 PT-1000-B-L-R	100 - 1000	Racked	No	96/rack, 10 racks/pack, 5 packs/case
78692-478 PT-1000-B-R-S	78692-474 PT-1000-B-L-R-S	100 - 1000	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case

1000 µL Wide Bore



--	78692-492 PT-1005-WB-C-L	100 - 1000	Bulk	No	1,000/pack, 10 packs/case
--	78692-494 PT-1005-WB-C-L-R	100 - 1000	Racked	No	96/rack, 10 racks/pack, 5 packs/case
--	78692-496 PT-1005-WB-L-R-S	100 - 1000	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case

1200 µL Fine point



78692-498 PT-1200-C	78692-500 PT-1200-C-L	1 - 1200	Bulk	No	1,000/pack, 10 packs/case
78692-506 PT-1200-R	78692-502 PT-1200-L-R	1 - 1200	Racked	No	96/rack, 8 racks/pack, 4 packs/case
78692-508 PT-1200-R-S	78692-504 PT-1200-L-R-S	1 - 1200	Racked	Yes	96/rack, 8 racks/pack, 4 packs/case

Empty Racks for Reloads

VWR Corning Cat. No.	Description	Sterile	Quantity
78692-458 PRFL-R-300	Empty Rack for 10 µL	No	10/pack, 5/case
78692-456 PRFL-R-200	Empty Rack for 200 µL	No	10/pack, 5/case
78692-454 PRFL-R-1000	Empty Rack for 1000 µL	No	10/pack, 5/case

Axygen® Macrovolume Pipette Tips

VWR Corning Cat. No.	Description	Pack Type	Sterile	Quantity
89005-638 T-5000-C	5 mL, clear	Bulk	No	250 tips/pack, 10 packs/case
89005-641 T-10ML-C	10 mL, clear	Bulk	No	100 tips/pack, 20 packs/case

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

For additional product or technical information, visit vwr.com/corning, call 1.800.932.5000, or contact your VWR representative.



vwr.com

Any prices, product, and/or services details are current when published and subject to change. Offers subject to product availability. Void where prohibited by law or company policy. Any claims concerning products manufactured by third parties are the sole claims of the manufacturer and not those of Avantor or its affiliates. All products are subject to the terms and conditions of sale of Avantor or its affiliates, including the applicable limitations of liability contained therein. The information herein should not be considered a warranty of any kind. Please carefully review any product information sheets and safety information for the product. Trademarks are owned by Avantor Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy and additional disclaimers. © 2026 Avantor Inc. All rights reserved.