

BRAND® HandyStep® Repeating Pipettes



BRAND HandyStep touch/touch S

Choose the HandyStep touch for your repetitive pipetting tasks—dispensing, auto-dispensing and the ability to store favorites. The HandyStep touch S adds such advanced features like multi-aspiration, sequential dispensing and titration. Paired with the PD-Tip™ // Dispenser Tips with tip size coding, they make the perfect team for easy pipetting.

- **Intuitive touch screen navigation:** Easily move through modes and settings. Plus integrated help function is right where it's needed.
- **Recognizes Tips Automatically:** The pipette automatically recognizes the size of encoded BRAND PD-Tip II Precision Dispenser Tips upon tip insertion and displays the tip's volume, default dispensing volume, and default number of dispensing steps.
- **Reduces Repetitive Motion by up to 97%:** The instrument features a unique Auto-Dispense operating mode in which the pipette “learns” the operator's dispensing rhythm and dispenses at that rhythm as long as the control button is depressed – no programming needed!

Use the tips you want

The only touch screen repeating pipette compatible with all standard syringe tips, the HandyStep touch accepts:

- BRAND PD-Tip II Precision Dispenser Tips
- Rainin® Encode™ syringe tips
- Eppendorf® Combitips advanced®
- Fisherbrand® Dispenser Tips
- VWR® Combi-Syringes
- Gilson® Repet-Tips

BRAND HandyStep S

The HandyStep S repeating pipette provides unprecedented versatility in a purely mechanical repeating pipette. In conjunction with BRAND PD-Tip II Precision Dispenser Tips, the HandyStep S offers nearly sixty different volume settings with accuracy, precision and dependability.

- **Comfortable:** The slim design of the HandyStep S is lightweight, well-balanced, and provides a convenient control layout for smooth dispensing and one-handed volume selection.
- **Reliable:** Robust construction, inspired design and innovative polymers deliver a dependable instrument.
- **Ready-To-Go:** A purely mechanical design with no batteries to charge or replace means the HandyStep S is ready when you are. A convenient integrated volume table provides a quick reference for dispensed volumes.
- **Flexible:** Ten sizes of BRAND PD-Tip II Precision Dispenser Tips allow the HandyStep S to dispense 59 different volumes. Open design allows use of third party tips—no more being locked into a sole source!
- **Contact-Free:** A tip ejection system allows the user to remove tips without handling.

	BRANDTECH	VWR
Description	Cat. No.	Cat. No.
HandyStep S	705110	89176-652
HandyStep touch	705200	76338-254
HandyStep touch S	705210	76338-256
HandyStep touch including charging stand	705201	76369-710
HandyStep touch S including charging stand	705211	76369-712
Charging stand, all Handystep touch models	705220	76338-258

BRAND® PD-Tip™ II Syringe Tips

For Repeating Pipettes

The Perfect Tips for Your Repeating Pipette

BRAND encoded PD-Tip II syringe tips are manufactured from high-quality virgin plastics for use with all standard repeating pipettes in a wide variety of applications.

- **Compatible with Most Pipettes:** BRAND PD-Tip II Precision Dispenser Tips are compatible with repeating pipettes using standard tips. BRAND PD-Tip II Precision Dispenser Tips are encoded for automatic recognition by the BRAND HandyStep touch repeating pipette, the Rainin® AutoRep™ E pipette, and the Gilson® Repetman®.
- **Reliable:** A precision seal between the piston and cylinder ensure smooth operation, increased accuracy, and increased precision. Extended piston design minimizes dead volume and maximizes fluid utilization.
- **Manufactured from High Quality Materials:** PD-Tip II syringe tips are made entirely from high-quality virgin polyethylene and polypropylene. (Piston of 0.1 mL tip is liquid crystal polymer-LCP).
- **BIO-CERT® PD-Tip II Precision Dispenser Tips** are sterile, individually wrapped, and certified free of DNA, RNase, endotoxins, and ATP to meet the most stringent sample integrity requirements.
- **Batch Certificate of Performance** which documents tip accuracy and precision testing is supplied with each pack of BRAND PD-Tip II Precision Dispenser Tips.



Capacity mL	Pack of	BRANDTECH non-sterile	VWR non-sterile	BRANDTECH sterile	VWR sterile
0.1	100	705700	76338-268	705730	76338-288
0.5	100	705702	76338-270	705732	76338-290
1.0	100	705704	76338-272	705734	76338-292
1.25	100	705706	76338-274	705736	76338-294
2.5	100	705708	76338-276	705738	76338-296
5	100	705710	76338-278	705740	76338-298
10	100	705712	76338-280	705742	76338-300
12.5	100	705714	76338-282	705744	76338-302
25	50 + 1 adapter	705716	76338-284	705746*	76338-304
50	25 + 1 adapter	705718	76338-286	705748	76338-306
Assorted	20/size	705720	76338-308	—	—
Adapter for PD-Tips II 25 / 50 mL tips	5	702398	18900-380	702399	47751-618

*BIO-CERT : Pack of 25 + 1 adapter




avantor™
delivered by **vwr™**

VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2024 Avantor, Inc. All rights reserved.