

Data sheet

Product: minicaps® 5 µl

Code-no.: 900 01 05

The above named product conforms with the requirements of the standard ISO 7550 and the IVD Directive 98/79/EG

Glass material: DURAN, borosilicate glass 3.3 (acc. ISO 3585)

Inner diameter: $0,45 \pm 0,02$ mm

Outer diameter: $0,89 \pm 0,02$ mm

Length: 31 ± 1 mm

Tolerances: Accuracy A: $\pm 0,5\%$
Coefficient of variation: $1,0\%$

Expiry date: 2,5 years after production date



.....
Benz, Quality manager/ 01.06.2021