



Certificate of Conformity

VWR Adhesive PCR Film Seal, Sterile

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No: 391-1255

Lot No: TG309F

Material: Polypropylene

Each film : L x D: 135.1 x 79.4 mm

Pack: 100

Sterility:

Product has been subjected to gamma radiation level appropriate to achieve 10^{-3} SAL.
Required dosage: 16.4kGy, Minimum dosage: 16.4kGy, Maximum dosage: 25kGy

Sterility date: 25th Nov 2019

Date of Manufacture: 05th Nov 2019

Product shelf life:

Storage recommendation: 3 years if maintained at 70 °F $\pm$ 10°; 50 % RH $\pm$ 10 %

Best if used by date: Nov 2022



Certificate of Conformity

We certify that the above stated product meets all specifications, dimensions, and tolerances as detailed in the product drawing, QC Visual & Performance specifications in QC procedures.

None of the substances listed below are knowingly or intentionally added during the manufacturing process and are not at or above 0.1% concentration level. The information presented was prepared based on the best knowledge and understanding of the Quality Control Dept.

Substance Name	Ec No.
Cadmium chloride	233-296-7
Diisohexyl phthalate	271-093-5
Sodium peroxometaborate	231-556-4
Sodium perborate	239-172-9 and 234-390-0

This product does not contain the following material as an intentionally added component.

Substance Name
Natural Rubber latex
Dry Natural Rubber
Animal Derivatives

This product has been passed final inspection and has been released for distribution.

This document has been produced electronically and is valid without a signature

Date of issue: **September 2020**