

Certificate of Conformity

VWR Catalogue Number	See List
Description	Centrifuge Tube with Flat Cap
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2025-07-08

Catalogue Number	Capacity [mL]	Colour	Sterility
525-1536	15	-	Nonsterile
525-1537	15	-	Nonsterile
525-1538	15	-	Sterile
525-1539	15	-	Sterile
525-1540	15	-	Sterile
525-1541	15	Amber	Nonsterile
525-1542	15	Amber	Nonsterile
525-1543	15	Amber	Sterile
525-1544	15	Amber	Sterile
525-1545	15	Amber	Sterile
525-1546	50	-	Nonsterile
525-1547	50	-	Nonsterile
525-1548	50	-	Sterile
525-1549	50	-	Sterile
525-1550	50	-	Sterile
525-1551	50	Amber	Nonsterile
525-1552	50	Amber	Nonsterile
525-1553	50	Amber	Sterile
525-1554	50	Amber	Sterile
525-1555	50	Amber	Sterile
525-1556	50	-	Sterile
525-1557	50	-	Sterile

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 & ISO 13485:2016**

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Material	Tube: Virgin Polypropylene (PP).
	Cap: Polyethylene (PE).
	Resins are USP Class-VI certified, FDA regulation CFR 21, and are FDA approved for food contact.
Temperature Resistance	-80 to +121 °C.
Sterilization	Product labelled as sterile are irradiated and dose released upon ISO 11137 recommended practices in effect at the time of validation meet a minimum requirement of 10 ⁻⁶ SAL (which is not applicable for Non-Sterile products).
ATP Assay	Not applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Endotoxin levels are tested by Limulus Amoebocyte Lysate (LAL) Assay. Product is free from Endotoxin.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	This product is free of Bisphenol A.
DEHP Statement	Please see RoHS information below.

RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Room Temperature.
Shelf Life	3 years from the date of manufacture (which is not applicable for Non-Sterile products).

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	04-2021
2	2.0	Modified Quality system compliance and disclaimer	07-2025