

Date: 26 January 2018
Quality Assurance VWR Chemicals

VWR International BVBA
Researchpark Haasrode 2020
Geldenaaksebaan 464
B - 3001 Leuven
e-mail: qaqc@vwr.com

Subject: VWRC24360.233 GELATINE POWDER PH.EUR. manufacturing flow

Dear Customer,

We hereby confirm that the manufacturing flow for VWRC24360.233 GELATINE POWDER PH.EUR is as follows:

Bovine bones are subjected to an acid treatment. The resulting ossein is treated with lime solution, washed and extracted. After purification via filtration and deionisation the resulting solution is concentrated and sterilized, the target product is extruded, dried, milled and packed.

In case of any remarks or questions, please do not hesitate to contact us.

Best regards,

QA Private Label Chemicals
Matthias Pues

