

# Certificate of Conformity

|                                   |                                        |
|-----------------------------------|----------------------------------------|
| <b>VWR Catalogue Number</b>       | See List                               |
| <b>Description</b>                | VWR®, Test tubes with rim, thin walled |
| <b>Country of Origin</b>          | Manufactured in India                  |
| <b>Date of Issue (yyyy-mm-dd)</b> | 2023-03-31                             |

| Catalogue Number | Length<br>[mm] | External Diameter<br>[mm] | Wall Thickness<br>[mm] |
|------------------|----------------|---------------------------|------------------------|
| 212-0753         | 70             | 8                         | 0.4 – 0.5              |
| 212-0759         | 75             | 10                        | 0.4 – 0.5              |
| 212-0760         | 75             | 12                        | 0.4 – 0.5              |
| 212-0754         | 100            | 10                        | 0.4 – 0.5              |
| 212-0755         | 100            | 12                        | 0.4 – 0.5              |
| 212-0761         | 130            | 14                        | 0.4 – 0.5              |
| 212-0762         | 130            | 16                        | 0.5 – 0.6              |
| 212-0763         | 150            | 25                        | 0.6 – 0.7              |
| 212-0756         | 160            | 16                        | 0.5 – 0.6              |
| 212-0757         | 180            | 18                        | 0.5 – 0.6              |
| 212-0758         | 180            | 20                        | 0.5 – 0.6              |
| 212-0764         | 200            | 25                        | 0.6 – 0.7              |
| 212-0765         | 200            | 30                        | 0.7 – 0.8              |

## Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard and calibrated in **ISO 17025:2017** certified Laboratory. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

## Product Specifications

The product is in conformance with the requirements of the following standards and regulations:

Technical standards for the products: NA

Calibration procedure: NA

|                                |                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Material</b>                | The products are made from glass tubing that conforms to ASTM E-438 Type-1, Class B Borosilicate Glass 5.1 Expansion.                                                                                                                                                                                                                                       |
| <b>Sterilization</b>           | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>ATP Assay</b>               | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>DNase &amp; RNase Free</b>  | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE). |
| <b>Cytotoxicity</b>            | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Non-Pyrogenic Statement</b> | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Latex Statement</b>         | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>BPA Statement</b>           | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>DEHP Statement</b>          | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>RoHS</b>                    | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>REACH Statement</b>         | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage Conditions</b>      | Room temperature.                                                                                                                                                                                                                                                                                                                                           |
| <b>Shelf Life</b>              | No predefined shelf life.                                                                                                                                                                                                                                                                                                                                   |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.

## Version History

| Sl.no | Version number | Description                          | Month & Year |
|-------|----------------|--------------------------------------|--------------|
| 1     | 1.0            | Created General CoQ using new layout | 03-2023      |
|       |                |                                      |              |
|       |                |                                      |              |