

Certificate of Conformity

VWR Catalogue Number	See List
Description	Prime overhead stirrers, VWR®
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2022-06-09

Catalogue Number	Description	Max. stirring capacity H ₂ O [L]	Max. Viscosity [mPas]	Speed range [min ⁻¹]
441-0510	VWR® Prime overhead stirrer 8 Ncm	8	4000	50 - 2000
441-0511	VWR® Prime overhead stirrer 15 Ncm	15	8000	50 - 1000
441-0508	VWR® Prime overhead stirrer 20 Ncm	20	10000	50 - 2000
441-0507	VWR® Prime overhead stirrer 40 Ncm	40	30000	50 - 1000

Quality System Compliance

Products are manufactured under the **ISO 9001 & ISO 13485** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Instrument Body Material - Polypropylene and Impeller - SS316L.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	Not Applicable.

Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Permissible Ambient Temperature: 5-40° C. Permissible Ambient Relative Humidity: 80%
Shelf Life	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.