

Material	20413.241
Material description	Nitric acid 68 %
Grade	GPR RECTAPUR
CAS Number	7697-37-2
Molecular formula	HNO ₃

Characteristics	Specifications
Assay	68.0 --> 70.0 %
Density (20/4)	1.400 --> 1.420
Heavy metals (as Pb)	Max. 20 ppm
Ignition residue (SO ₄)	Max. 50 ppm
Cl (Chloride)	Max. 5 ppm
SO ₄ (Sulphate)	Max. 5 ppm
Fe (Iron)	Max. 2 ppm
Pb (Lead)	Max. 1 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250 Briare;
France