

|                      |                            |
|----------------------|----------------------------|
| Material             | 23857.320                  |
| Material description | Petroleum spirit 100-140°C |
| Grade                | GPR RECTAPUR               |

|            |            |
|------------|------------|
| CAS Number | 64742-49-0 |
|------------|------------|

| Characteristics                                       | Specifications |
|-------------------------------------------------------|----------------|
| Substances coloured by H <sub>2</sub> SO <sub>4</sub> | Passes test    |
| Acidity                                               | ≤ 0.0003 meq/g |
| Density (20/4)                                        | 0.690 - 0.750  |
| Distillation (min. 90 %)                              | 100 - 140 °C   |
| Evaporation residue                                   | ≤ 20 ppm       |

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare  
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250  
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.