

Certificate of Conformity

VWR Catalogue Number	See List
Description	Sampling spoons, Colour white
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2026-01-21

Catalogue Number	Capacity [mL]	Length [mm]	Width [mm]
231-0493	2.5	127	30
231-0494	10	170	41.5

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with [ISO 9001:2015](#).

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polystyrene (PS). Materials used to produce the products are listed in List of Authorised Substances of Regulation EU/10/2011. The material grades conform to EU Directive 1935/2004 on materials and articles intended to come into contact with food.
Sterilization	Products has been gamma irradiated and dose released upon ISO11137-1 which meets a minimum requirement of 10^{-6} SAL with a specified dose range of 17.5-50 kGy.
ATP Assay	Not Tested.
DNase & RNase Free	Not Tested.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Tested.
Non-Pyrogenic Statement	Not Tested.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Store in a room temperature.
Shelf Life	4 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	08-2021
2	1.1	Width included in table	11-2023
3	2.0	Modified Quality system compliance and disclaimer	07-2025
4	2.1	Shelf Life updated as 4 years	01-2026