

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Pasteur Pipettes
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2025-07-16

Catalogue Number	Description	Capacity [mL]	Length [mm]
612-1701	Non plugged	2	150
612-1702	Non plugged	2	230
612-3850	Non plugged, sealed tip, breakable	2	230
612-3860	Non plugged, sealed tip, breakable	2	270
612-1798	Pre-plugged	2	150
612-1799	Pre-plugged	2	230
612-3849	Pre-plugged, sealed tip, breakable	2	230
612-3861	Pre-plugged, sealed tip, breakable	2	270

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with [ISO 9001:2015](#).

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Top body outer diameter: 7.1 mm

Wall thickness: 0.53 mm

Jet outer diameter: 1.5 mm

Material	Soda Lime glass tubing.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	Not Applicable.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Storage at ambient temperature.
Shelf Life	Not Applicable.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	10-2023
2	2.0	Modified Quality system compliance and disclaimer	07-2025