




Dimensions (mm)		
Code		
D1	27.30	+/- 0.20
D5	21.65	+/- 0.25
D6	19.40	+/- 0.30
D7	16.00	+ 0.50
H1	57.50	+/- 0.50
H3	9.50	+/- 0.30
Thread Turns	Min - 1.50	
Neck Finish		
Neck Type	Screw Neck	
Closure Type	22mm R3	
Materials		
Tube	Clear 33	
Closure	Expansion Urea	
Packing		
No/Carton	500	

T006/01

Scintillation Vial

Low Potassium Glass



This drawing is property of S Murray & Co

Issue: G3

Date: April 2023

Signed: Geoff Hambly, Quality Director

