


Recombinant Human Thrombopoietin Protein (6xHis Tag) (A317619)

Specifications:

Name:	Recombinant Human Thrombopoietin Protein (6xHis Tag)
Applications:	SDS-PAGE
Expression System:	HEK293 cells
Nature:	Recombinant
Protein Species:	Human
Protein Length:	Protein fragment.
Amino Acids:	Ser22-Gly353
Tag:	N-terminal 6xHis Tag and C-terminal 6xHis Tag
Molecular Weight:	37.3 kDa
Conjugate:	Unconjugated
Purity:	> 95% by SDS-PAGE under reducing conditions.
Product Form:	Lyophilized
Concentration:	Reconstitution dependent.
Formulation:	Lyophilized from a 0.2 μ m filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
Storage:	Shipped at 4°C. Lyophilized: Store at -20°C to -80°C. Reconstituted: Aliquot and store at -80°C. Product is stable for one year. Avoid freeze/thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

Images:



Recombinant Human Thrombopoietin Protein (6xHis Tag) (A317619) has greater than 95% purity as determined by SDS-PAGE under reducing conditions.