




Tamper Evident Pull type Cap- 33mm

Tamper Evident Threaded Cap- 45mm

Tamper Evident Threaded Cap- 65mm



Bottle affixed with label

- Unique Barcode for easy traceability
- Sterility validity
- Lot #
- Sodium Thiosulphate Dosage
 - Blue = UNDOSED
 - RED = 20mg / L
 - GREEN = 120mg / L
- Write on area for sample information

CAUTION: DO NOT RINSE THE BOTTLE

Material

Bottle : HDPE, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.

Closure : Polypropylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed With Food grade colour.

BSE/TSE Statement: These products are manufactured without the use of Raw materials of animal origin.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.



Material Recycle Code



Description: Water Sampling Bottle (HDPE)

Drawing:
Not to Scale



Release: 2022/12/05

Revision: 1

Cat#	Product Description	Capacity (ml)	Sodium Thiosulphate Dosage	Material (Bottle / Cap)	Bottle Shape	Sterile	Colour (Bottle / Cap)	Neck Size	Neck (D)(mm)	Length (mm)	Width (mm)	Height without Cap (mm)	Height with Cap (mm)	General Wall Thickness (mm)	Pack Unit	Case Size
331-0383	HDPE BOTTLE 500ML UNDOSED 33mm Cap ST	500 ml	Undosed	HDPE / PP	Square	Yes	Natural / RED	33mm	29mm	71.5	71.5	149.5	150.5	0.7	CS	100 EACH
331-0384	HDPE BOTTLE 500ML 20mg/L 33mm Cap ST	500 ml	20mg/L	HDPE / PP	Square	Yes	Natural / RED	33mm	29mm	71.5	71.5	149.5	150.5	0.7	CS	100 EACH
331-0385	HDPE BOTTLE 500ML 120mg/L 33mm Cap ST	500 ml	120mg/L	HDPE / PP	Square	Yes	Natural / RED	33mm	29mm	71.5	71.5	149.5	150.5	0.7	CS	100 EACH
331-0386	HDPE BOTTLE 1000ML UNDOSED 33mm Cap ST	1000 ml	Undosed	HDPE / PP	Square	Yes	Natural / RED	33mm	29mm	83	83	197	199	0.7	CS	84 EACH
331-0387	HDPE BOTTLE 1000ML 20MG/L 33mm Cap ST	1000 ml	20mg/L	HDPE / PP	Square	Yes	Natural / RED	33mm	29mm	83	83	197	199	0.7	CS	84 EACH
331-0388	HDPE BOTTLE 1000ML 120MG/L 33mm Cap ST	1000 ml	120mg/L	HDPE / PP	Square	Yes	Natural / RED	33mm	29mm	83	83	197	199	0.7	CS	84 EACH
331-0389	HDPE BOTTLE 250ML UNDOSED 45mm cap ST	250 ml	Undosed	HDPE / PP	Square	Yes	Natural / RED	45mm	36mm	56.5	56.5	108	110	0.7	CS	210 EACH
331-0390	HDPE BOTTLE 250ML 20mg/L THIO 45mm Cap ST	250 ml	20mg/L	HDPE / PP	Square	Yes	Natural / RED	45mm	36mm	56.5	56.5	108	110	0.7	CS	210 EACH
331-0391	HDPE BOTTLE 500ML UNDOSED 45mm Cap ST	500 ml	Undosed	HDPE / PP	Square	Yes	Natural / RED	45mm	36mm	70	70	146.5	148.5	0.7	CS	100 EACH
331-0392	HDPE BOTTLE 500ML 20MG/L 45mm Cap ST	500 ml	20mg/L	HDPE / PP	Square	Yes	Natural / RED	45mm	36mm	70	70	146.5	148.5	0.7	CS	100 EACH
331-0393	HDPE BOTTLE 500ML 120MG/L 45mm ST	500 ml	120mg/L	HDPE / PP	Square	Yes	Natural / RED	45mm	36mm	70	70	146.5	148.5	0.7	CS	100 EACH
331-0394	HDPE BOTTLE 1000ML UNDOSED 45mm ST	1000 ml	Undosed	HDPE / PP	Square	Yes	Natural / RED	45mm	36mm	83	83	193	195	0.7	CS	84 EACH
331-0395	HDPE BOTTLE 1000ML 20MG/L 45mm ST	1000 ml	20mg/L	HDPE / PP	Square	Yes	Natural / RED	45mm	36mm	83	83	193	195	0.7	CS	84 EACH
331-0396	HDPE BOTTLE 1000ML 120MG/L 45mm ST	1000 ml	120mg/L	HDPE / PP	Square	Yes	Natural / RED	45mm	36mm	83	83	193	195	0.7	CS	84 EACH
331-0397	HDPE BOTTLE 500ML UNDOSED 63mm ST	500 ml	Undosed	HDPE / PP	Square	Yes	Natural / RED	63mm	56mm	75.5	75.5	125	127	0.7	CS	160 EACH
331-0398	HDPE BOTTLE 500ML 20MG/L 63mm ST	500 ml	20mg/L	HDPE / PP	Square	Yes	Natural / RED	63mm	56mm	75.5	75.5	125	127	0.7	CS	160 EACH
331-0399	HDPE BOTTLE 500ML 120MG/L 63mm ST	500 ml	120mg/L	HDPE / PP	Square	Yes	Natural / RED	63mm	56mm	75.5	75.5	125	127	0.7	CS	160 EACH
331-0400	HDPE BOTTLE 1000ML UNDOSED 63mm ST	1000 ml	Undosed	HDPE / PP	Rectangular	Yes	Natural / RED	63mm	56mm	94	74.5	182	184	0.7	CS	90 EACH
331-0401	HDPE BOTTLE 1000ML 20MG/L 63mm ST	1000 ml	20mg/L	HDPE / PP	Rectangular	Yes	Natural / RED	63mm	56mm	94	74.5	182	184	0.7	CS	90 EACH
331-0402	HDPE BOTTLE 1000ML 120MG/L 63mm ST	1000 ml	120mg/L	HDPE / PP	Rectangular	Yes	Natural / RED	63mm	56mm	94	74.5	182	184	0.7	CS	90 EACH
331-0403	HDPE BOTTLE 250ML 20MG/L THIO IND.PACK	250ml	20mg/L	HDPE / PP	Square	Yes	Natural / RED	45mm	36mm	56.5	56.5	108	110	0.7	CS	210 EACH
331-0404	HDPE BOTTLE 500ML 20MG/L THIO IND.PACK	500ml	20mg/L	HDPE / PP	Square	Yes	Natural / RED	45mm	36mm	70	70	146.5	148.5	0.7	CS	130 EACH
331-0405	HDPE BOTTLE 500ML 120MG/L THIO IND.PACK	500ml	120mg/L	HDPE / PP	Square	Yes	Natural / RED	45mm	36mm	70	70	146.5	148.5	0.7	CS	130 EACH
331-0406	HDPE BOTTLE 500ML 20MG/L THIO IND.PACK	500ml	20mg/L	HDPE / PP	Square	Yes	Natural / RED	45mm	36mm	70	70	146.5	148.5	0.7	CS	100 EACH
331-0407	HDPE BOTTLE 500ML 120MG/L THIO IND.PACK	500ml	120mg/L	HDPE / PP	Square	Yes	Natural / RED	45mm	36mm	70	70	146.5	148.5	0.7	CS	100 EACH
331-0408	HDPE BOTTLE 1L 20MG/L THIO IND.PACK	1000ml	20mg/L	HDPE / PP	Rectangular	Yes	Natural / RED	63mm	56mm	94	74.5	182	184	0.7	CS	77 EACH
331-0409	HDPE BOTTLE 1L 120MG/L THIO IND.PACK	1000ml	120mg/L	HDPE / PP	Rectangular	Yes	Natural / RED	63mm	56mm	94	74.5	182	184	0.7	CS	77 EACH