


Purification made simple

A complete workflow solution from sample collection to results


Automate your workflow in 4 simple steps





Automate at Any Scale

EXTRACTION MADE SIMPLE AND CONSISTENT



VWR® 96 SAMPLE PURIFICATION SYSTEM

- On-board touchscreen programming
- 8-position carousel with 2 heating positions
- Sample volume 50 µl to 1000 µl
- Process up to 96 samples at a time

Description	Cat. No.
VWR® 96 Sample purification system 115 V	77321-884
VWR® 96 Sample purification system 230 V	77321-886
VWR® 96 Tips combs CS50	77321-882
J.T.Baker® Deep well plate SRW magnetic	43001-0504



VWR® 16 SAMPLE PURIFICATION SYSTEM


- Compact design
- Process up to 16 samples at a time
- Programmable with PC software or Android app
- Sample volume 50 µl to 1000 µl

Description	Cat. No.
VWR® 16 Sample purification system 115 V	77321-876
VWR® 16 Sample purification system 230 V	77321-878
VWR® Tablet	77321-874
VWR® 8 Place tips combs CS50	77321-872
J.T.Baker® Deep well plate SRW magnetic	43001-0504



Explore the collection

COLLECTION SOLUTIONS DESIGNED FOR ANY SAMPLE



DNA/RNA SHIELD

- Stable transport at room temperature of DNA & RNA
- Compatible with cells, tissue, blood, fluids, FFPE, stool, soil, etc.
- Devices formatted for major automation platforms
- Safe handling & pathogen inactivation


Description	Cat. No.
DNA/RNA Shield™ (1X Concentrate)	76020-418
DNA/RNA Shield™ SafeCollect Swab Kit 2 ml	76628-454
DNA/RNA Shield™ SafeCollect Saliva Kit 2 ml	76628-464
DNA/RNA Shield™ Collection Tube with Swab (2 ml fill)	76020-432
DNA/RNA Shield™ Fecal Collection Tube (with beads)	77020-304
DNA/RNA Shield™ Blood Collection Tube	76020-434



Discover quality purification

AUTOMATED WORKFLOWS FOR NGS-READY NUCLEIC ACID PURIFICATION

High-quality DNA recovery




QUICK PURIFICATION KITS

- High-quality nucleic acid recovery from any sample (cells, tissue, blood, fluids, FFPE, stool, soil, environmental)
- NGS-ready nucleic acid ready for any application
- DNase I Included where applicable

Description	Cat. No.
Quick-DNA™ /RNA HT	77322-178
Quick-DNA™ Magbead Plus Kit	76197-100
Quick-DNA™ Magbead	76421-026

High-quality total RNA



DIRECT-ZOL-96 MAGBEAD RNA

- Bypass chloroform, phase separation and precipitation steps
- Ultra-pure RNA, no phenol or DNA carryover
- Complete RNA recovery without miRNA loss

Description	Cat. No.
Direct-zol™-96 MagBead RNA Kit	76211-322



See the Full Picture

THE ALL-IN-ONE PACKAGE FOR PURITY ANALYSIS



VWR® NANO SPECTROPHOTOMETER - THE ALL-IN-ONE PACKAGE FOR PURITY ANALYSIS

- Rapid results from a single droplet
- Extended-life Xenon lamp & silicon CCD array for full spectral scanning
- Built-in printer and cuvette port for compact, all-in-one footprint
- Detects up to 4,500 ng/μl of undiluted dsDNA

Description	Cat. No.
VWR® Nano Spectrophotometer, Drop and Cuvette	76777-150



Support from automation specialists
and ready-to-run scripts for all kits
available at no extra cost



Scan to Inquire About Ordering Information or
Additional Customer Support

Or Call 1-800-932-5000 | 7:30 a.m. to 8:00 p.m. (EST)