



Sol - Vet[®] Insulin Syringes

- Needles are tri-bevel, anti-coring and silicone oil lubricated for a smoother injection
- Syringes require less penetration force, resulting in greater animal comfort
- Low dead space results in over 90% less medication waste compared to the ISO standard for conventional insulin syringes
- Available in U-40 and U-100
- Bold, red unit markings to indicate U40 insulin syringe with half unit options available for accurate dosing, U100 is indicated by orange cap
- Large flange offers secure grip and easier medication dispensing

***samples available upon request**

VWR Cat No.

Product Description

76582-940	Sol-Vet 1mL U-40 Insulin Syringe w/ Fixed Needle 29G x 1/2" (PE Bag)
76582-944	Sol-Vet 1mL U-40 Insulin Syringe w/ Fixed Needle 30G x 1/2" (PE Bag)
76582-956	Sol-Vet 0.5mL U-40 Insulin Syringe w/ Fixed Needle 29G x 1/2" (half unit markings)(PE Bag)
76582-960	Sol-Vet 0.5mL U-40 Insulin Syringe w/ Fixed Needle 30G x 1/2" (half unit markings)(PE Bag)
76582-946	Sol-Vet 0.3mL U-40 Insulin Syringe w/ Fixed Needle 29G x 1/2 (half unit markings)(PE Bag)
76582-952	Sol-Vet 0.3ml U-40 Insulin Syringe w/ Fixed Needle 30G x 1/2 (half unit markings)(PE Bag)
76582-938	Sol-Vet 1mL U-100 Insulin Syringe w/ Fixed Needle 29G x 1/2" (PE Bag)
76582-942	Sol-Vet 1mL U-100 Insulin Syringe w/ Fixed Needle 30G x 1/2" (PE Bag)
76582-954	Sol-Vet 0.5mL U-100 Insulin Syringe w/ Fixed Needle 29G x 1/2" (half unit markings)(PE Bag)
76582-958	Sol-Vet 0.5mL U-100 Insulin Syringe w/ Fixed Needle 30G x 1/2" (half unit markings)(PE Bag)
76582-948	Sol-Vet 0.3mL U-100 Insulin Syringe w/ Fixed Needle 29G x 1/2" (half unit markings)(PE Bag)
76582-950	Sol-Vet 0.3mL U-100 Insulin Syringe w/ Fixed Needle 30G x 1/2" (half unit markings)(PE Bag)