




PRODUCT IMAGE FOR
125ml, 250ml, 500ml & 1000ml
CAPACITY BOTTLES



PRODUCT IMAGE FOR
2000ml CAPACITY BOTTLE

REVISION NO	REVISION DATE	DESCRIPTION	REVISED BY	APPROVED BY

RECTANGULAR WIDE MOUTH BOTTLE SPECIFICATION

VWR-US CAT NO	VWR EUROPE CAT NO	MATERIAL	COLOR	CAPACITY (ml)	APROX BRIM CAPACITY (ml)	CAP SIZE	HEIGHT WITH OUT CAP (mm)	BODY WIDTH (mm)	BODY LENGTH (mm)
414004-191	215-1905	HDPE	TRANSLUCENT	125ml	140ml	28-415	97	38	59
414004-192	215-1906	HDPE	TRANSLUCENT	250ml	286ml	38-415	112	51.5	74.5
414004-193	215-1907	HDPE	TRANSLUCENT	500ml	620ml	48-415	114	64	99.5
414004-194	215-1908	HDPE	TRANSLUCENT	1000ml	1130ml	53-415	180	67	127
414004-195	215-1909	HDPE	TRANSLUCENT	2000ml	2340ml	63	237	96	150
414004-196	-	HDPE	AMBER	125ml	140ml	28-415	97	38	59
414004-197	-	HDPE	AMBER	250ml	286ml	38-415	112	51.5	74.5
414004-198	-	HDPE	AMBER	500ml	620ml	48-415	114	64	99.5
414004-199	-	HDPE	AMBER	1000ml	1130ml	53-415	180	67	127
414004-200	-	HDPE	AMBER	2000ml	2340ml	63	237	96	150

NOTE : 2L PRODUCT HAS INTEGRATED HANDLE DESIGN.

Closure Application Torque		
Closure Size	Minimum Torque Kg cm	Maximum Torque Kg cm
13-415	6	8
20-415	11	16
24-415	13	19
28-415	18	25
38-415	25	35
38-430	31	38
43-415	32	44
48-415	34	48
53-415	38	52
53T	43	60
63-415	46	64
83T	69	96
100-415	54	74

TITLE : WIDE MOUTH RECTANGULAR BOTTLE

DRAWN BY	
DRAWN DATE	2/2/2009
DRAWING NO	

Material : AS PER TABLE

Rev : 0 SCALE : NOT TO SCALE

