



Certificate of Conformity

VWR DISPOSABLE BEAKER CUPS

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

VWR Cat. No.	Height mm	Bottom diameter mm	Top inner diameter mm	Top outer diameter mm	Capacity ml	Pk
611-1363	23	21	25	33	5	1.000
611-1364	30	22	29	35	10	1.000
611-1365	32	29	43	49	20	500
611-1366	50	34	53	60	50	500
611-1367	65	56	69	82	150	100
611-1368	88	66	80	91	250	100
611-1369	124	101	121	148	800	100
611-1370	118	113	123	135	1000	100

Material : Polystyrene

Other Characteristics:

- Economical, multipurpose beaker for a very wide range of uses.
- Suitable for use in various sample dispensers
- Beakers with capacity of 5 and 20 ml are recommended for use in flame photometric analyses
- Excellent chemical resistance to dilute or weak acids, aqueous solutions, alcohols and bases.
- Excellent for use with atomic absorption, flame photometers and blood cell counters.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **November 2015**