



Specification

1.05978.0500 OGYE agar, Base acc. to ISO for microbiology

Specification	
Appearance (clearness)	clear
Appearance (colour)	yellow
pH-value (25 °C)	6.4 - 6.8

Specification		
Recovery on test medium (Candida albicans ATCC 10231 (WDCM 00054))	≥ 50	%
Recovery on test medium (Saccharomyces cerevisiae ATCC 9763 (WDCM 00058))	≥ 50	%
Recovery on test medium (Aspergillus brasiliensis (formerly A. niger) ATCC 16404 (WDCM 00053))	≥ 50	%
Recovery on test medium (Penicillium aurantiogriseum ATCC 16025)	≥ 50	%
Growth (Escherichia coli ATCC 25922 (WDCM 00013))	none	
Growth (Bacillus subtilis ATCC 6633 (WDCM 00003))	none	

Stefanie Fischer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.