



## VWR® CHEMICALLY RESISTANT APRONS

As per the manufacturer, the below product meets the following criteria:

|                             |                             |
|-----------------------------|-----------------------------|
| <b>VWR Catalog Number:</b>  | 89232-928                   |
| <b>Lot Number:</b>          | 10191915                    |
| <b>Description:</b>         | VWR BASIC APRON YELLOW CS50 |
| <b>Date of Manufacture:</b> | December 13, 2019           |
| <b>Country of Origin:</b>   | CHINA                       |

### Quality System Compliance

Products are manufactured under the ISO 9001:2015 & ISO13485:2003 standard. Products are Inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP. Inspection records are reviewed and signed off by qualified personnel for product release.

### QC Testing

Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E

|                                        |                           |            |            |
|----------------------------------------|---------------------------|------------|------------|
| <b>-Break Tensile properties(N):</b>   | ASTM D5034standard,       | MD= 181.20 | CD= 167.30 |
| <b>-Elongation rate (%):</b>           | ASTM D5034standard,       | MD= 62.50  | CD= 56.20  |
| <b>-Trapezoid tearing strength(N):</b> | FZ/T 01085-1999 standard, | MD= 816.28 | CD= 769.34 |

### Product Specifications

|                                |                                                                                                                                                                                                                                               |
|--------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Material:</b>               | SPP+PE 88g/m <sup>2</sup> Yellow (Spunbonded Polypropylene fabric+polyethylene fabric)                                                                                                                                                        |
| <b>BSE/TSE Statement:</b>      | The above product contains no ingredients of animal origin and no material derived from, or exposed to animals affected by or under quarantine for Transmitting, Animal Spongiform Encephalopathy (TSE)/Bovine Spongiform Encephalopathy(BSE) |
| <b>Non-Pyrogenic:</b>          | NA                                                                                                                                                                                                                                            |
| <b>DNase &amp; RNase Free:</b> | NA                                                                                                                                                                                                                                            |
| <b>Cytotoxicity:</b>           | NA                                                                                                                                                                                                                                            |
| <b>Latex Free Statement:</b>   | The above product is latex free.                                                                                                                                                                                                              |

Signed:



Jamie Ethier  
VP Global Quality  
VWR, Part of Avantor

**Date:** January 04, 2020