



Saniclave-Autoclave 50 - NON-FDA



Area of Focus Microbiology Veterinary

General lab use Academia Applications: Bacteria slides, Liquid media, Used Petri plates

Specifications

Pressure: 15 psi

Internal Temperature: 121°C

Sterilization Cycle: 30 minutes

Capacity: 8L

Net Weight: 5.4kg (12lbs)

Power Requirements: 110-130 VAC, 50-60 Hz, 1250 W

Display: Digital time and/or temperature

Durability: Rugged stainless steel chamber and cover

Safety: Over temperature alarm

Load type: Solid, non-lumened, unwrapped medical instruments

Features One button operation; Durable construction; Safe burn free exterior;

Weight *18 lbs*

Dimensions *17 × 17 × 17 in*