



## Product Specification Sheet

# VWR LABOTIMER DIGIT 100 MINUTES SMALL

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

*VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.*

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Cat. No: **609-0184**

Size: **46×89×12.7 mm**

Weight: **35.4 g**

Other Specifications:

- Easy-to-use timer has a repeat memory function for repetitive tests.
- Three-key operation can be learned in less than 10 seconds.
- May be placed on incubators, refrigerators, and ovens with the magnet on the back.
- Counts down from 99 minutes, 59 seconds in 1-second increments.
- Memory retains last programmed time, eliminating the need to re-enter test times.
- Alarm may be silenced manually or turns off automatically after 30 seconds.
- Supplied with silver-oxide battery, magnet, and Velcro®.

**Disclaimer:** VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **December 2012**

VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>



Version: 2.0