

Certificate of Conformity

VWR Catalogue Number	See List
Description	Erlenmeyer flasks, with standard ground joint
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2022-01-26

Flasks

Catalogue Number	Capacity [mL]	Height [mm]	Ext.Diameter [mm]	Socket Size [NS]
201-1385	25	70	42	14/23
201-1386	50	85	51	14/23
201-1362	50	85	51	24/29
201-1387	50	85	51	29/32
201-1388	100	105	84	14/23
201-1389	100	105	84	19/26
201-1363	100	105	84	24/29
201-1390	100	105	84	29/32
201-1369	200	131	79	24/29
201-1391	200	131	79	29/32
201-1370	250	140	85	24/29
201-1392	250	140	85	29/32
201-1372	300	156	87	24/29
201-1393	300	156	87	29/32
201-1373	500	175	105	24/29
201-1394	500	175	105	29/32
201-1395	1000	220	131	29/32

Hollow glass stoppers

Catalogue Number	NS
217-0887	10/19
217-0888	12/21
217-0889	14/23
217-0890	19/26
217-0891	24/29
217-0892	29/32
217-0893	34/35
217-0894	45/40

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard and calibrated in **ISO 17025:2017** certified Laboratory. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

The product is in conformance with the requirements of the following standards and regulations

Technical standards for products: Flask: DIN 12353.
Stopper: DIN 12252.

Calibration procedure: NA

Material Borosilicate 3.3 Glass DIN ISO 3585:1998 TYPE-1, CLASS A, and meets ASTM E-438, TYPE-1, CLASS-A Standard.

Sterilization Hot air sterilization, in the oven up to 140 °C
Steam sterilization, in an autoclave 121 °C/ 20 min
134 °C/ 10 min.

ATP Assay Not Applicable.

DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Store at room temperature in normal warehouse conditions.
Shelf Life	No predefined shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.