



Saving and storing, Reagent bottles

Reagent bottles, PP, GL 45



Transparent.

With screw cap made of PP.

Good chemical resistance, ideal for storage of liquids.

Autoclavable at 121 °C (2 bar) in compliance with DIN EN 285, except for size 5000 ml.

Suitable for contact with foodstuffs according to regulation (EU) No. 10/2011.

Cat. No. 100889: with handle, PE-HD.

ASK A QUESTION

KNOWLEDGE: SUITABILITY FOR FOODSTUFFS



Types / Sizes

Volume ml	Thread	Ø mm	Height mm	PU	Cat. No.
250	GL 45	73	132	20	101689
500	GL 45	87	172	10	101789
1000	GL 45	105	197	10	102089
2000	GL 45	131	241	6	102189
5000	GL 45	178	315	1	100889