




S:\DWGS\ISW900200

THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION AND IS A COPYRIGHT OF DWK LIFE SCIENCES INC. IT MAY NOT BE DUPLICATED FOR ANY PURPOSE IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN CONSENT OF DWK LIFE SCIENCES INC. ALL RIGHTS RESERVED. ELECTRONIC OR PAPER COPIES OF THIS PRINT OUTSIDE OF WHEATON'S SYSTEM ARE UNCONTROLLED AND ARE FOR REFERENCE ONLY UNLESS OTHERWISE STAMPED.

THIRD ANGLE PROJECTION



ISOMETRIC VIEW
SCALE 1:2



MATERIALS: DISH, CLEAR TYPE III SODA-LIME GLASS.
RACK, CLEAR TYPE III SODA-LIME GLASS.
HANDLE, #303-304 STAINLESS STEEL, OD 2mm.
COVER, CLEAR TYPE III SODA-LIME GLASS.

2. APPROX. WEIGHTS: DISH, 0.38kg; RACK, 0.15kg; HANDLE, 0.27oz; COVER, 0.16kg.

1. DIMENSIONS ARE IN MILLIMETERS.

---NOTES---

PRODUCT COMPONENTS		
ITEM	QTY	DESCRIPTION
1.	6	DISH
2.	6	RACK
3.	6	HANDLE
4.	6	COVER

TITLE 20 SLIDE STAINING DISH COMPLETE		DATE DRAWN 5/27/93	PROJ -	REVISED TO DWK STDS 09/16/19 EL
DRAWN BY BILL P. SHIPLEY		APPROX WEIGHT SEE NOTE #3	J	2016060815110173 9/27/16 TS
CAPACITY 20 SLIDE UNIT		CHECKED BY BCS	AUTHORIZED RB	REV ECN # DATE
FINISH		SIZE		C
SCALE: FULL UNLESS NOTED. DO NOT SCALE DRAWING.		DURAN WHEATON KIMBLE 1501 North 10th Street Milville, NJ 08332-2038 USA		