



designed for scientists

Technical Data



Technical Data

Type of movement	rolling
Permissible ON time	100 %
Speed min.	0 rpm
Speed fix	30 rpm
Speed display	none
Timer display	none
Operating mode	continuous operation
Rolls, number	10
Rolls, Ø	32 mm
Rolls, length	327 mm
Working area length	350 mm
Rolls, luffing angle fixed	3 °
Rolls, max stroke at roll end	16 mm
Rolls, changeable, removable	yes
Dimensions (W x H x D)	380 x 115 x 545 mm
Weight	8.242 kg
Permissible ambient temperature	4 - 50 °C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 21
Voltage	100 - 240 V
Frequency	50/60 Hz
Power input	24 W

DC Voltage	24 V=
Current consumption	1000 mA

IKA Shakers

ROLLER 10 basic

PRICE

INR 1,40,564

Ident. No.: 0004012000


1

V

ADD TO SHOPPING CART



< PRODUCT OVERVIEW



IKA Shakers

/// COMPARE MODELS NOW

COMPARE

