



VWR Maximum Protection Antifog Full Face Shields

- Protect Front and Sides of Face from Fluid Splash
- Low-Glare, Optically Clear PET Shield Allows Distortion-Free Vision

VWR Maximum Protection Antifog Full Face Shields provide total facial splash protection and can be worn for extended periods. Manufactured from PET (polyethylene terephthalate), these face shields feature an antifog coating and fit easily over eyewear. A soft, polyurethane foam sponge fits snugly against the forehead. Latex-free. Packaged 10 face shields per pack, five packs per box, four boxes per case. Dimensions: 33.5L x 21Hcm (13³/₁₆ x 8⁵/₈").

Face shields are manufactured in an ISO Certified facility under stringent process controls to ensure that each product meets exacting quality standards.

Face Shield Dimensions

Height (A)	Length (B)	Sponge Height (C)
8 ⁵ / ₈ "	13 ³ / ₁₆ "	1"

Technical Data

Test Item	Result	Test Standard
Thickness (mm)	0.18±0.02	ASTM D1777

VWR Maximum Protection Antifog Full Face Shields

Head Band	Cat. No.	Case of
Elastic	414004-679	200
Velcro	414004-866	200

To order, contact your VWR Sales Representative, visit vwr.com, or call **1.800.932.5000** today.