

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® LTS Compatible Pipet Tips
Country of Origin	Manufactured in Mexico
Date of Issue (yyyy-mm-dd)	2021-04-29

Racked, Filtered Tips

Catalogue Number	Volume [µL]	Length in mm (inches)	Sterility
76322-158	0.1 - 20	49.63 (1.95")	Sterile
76322-160	20 - 200	50.75 (2.00")	Sterile
76322-162	20 - 300	58.93 (2.32")	Sterile
76322-164	100 - 1000	88.9 (3.50")	Sterile
76327-216	100 - 1200	97.66 (3.85")	Sterile

Racked, non-filtered Tips

Catalogue Number	Volume [µL]	Length in mm (inches)	Sterility
76323-932	20 - 300	58.93 (2.32")	Sterile
76323-934	100 - 1000	88.9 (3.50")	Sterile
76323-944	0.1 - 20	49.63 (1.95")	Sterile
76323-946	20 - 200	50.75 (2.00")	Sterile
76327-218	100 - 1200	97.66 (3.85")	Sterile
76323-936	0.1 - 20	49.63 (1.95")	Non - sterile
76323-938	20 - 200	50.75 (2.00")	Non - sterile
76323-940	20 - 300	58.93 (2.32")	Non - sterile
76323-942	100 - 1000	88.9 (3.50")	Non - sterile
76327-220	100 - 1200	97.66 (3.85")	Non - sterile

Reload, non-filtered Tips

Catalogue Number	Volume [µL]	Length in mm (inches)	Sterility
76322-180	100 - 1000	88.9 (3.50")	Sterile
76322-530	0.1 - 20	49.63 (1.95")	Sterile
76322-532	20 - 200	50.75 (2.00")	Sterile
76322-534	20 - 300	58.93 (2.32")	Sterile
76327-212	100 - 1200	97.66 (3.85")	Sterile
76322-172	0.1 - 20	49.63 (1.95")	Non - sterile
76322-174	20 - 200	50.75 (2.00")	Non - sterile
76322-176	20 - 300	58.93 (2.32")	Non - sterile
76322-178	100 - 1000	88.9 (3.50")	Non - sterile
76327-210	100 - 1200	97.66 (3.85")	Non - sterile

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polypropylene.
Sterilization	Product labelled as sterile are irradiated and dose released upon ISO 11137-2:2013 recommended practices in effect at the time of validation meet a minimum requirement of 10 ⁻³ SAL (which is not applicable for non-sterile products).
ATP Assay	Products are tested for the presence of ATP and the assay resulted with a value of 7.3 ⁻⁹ mg/µL.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).

Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	No information available.
Latex Statement	The product is Latex free.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Store at ambient condition.
Shelf Life	See product label.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.