



CERTIFICATE OF ANALYSIS

LopHene® II

Lot # 0330201

Date of Manufacture March 2020

Expiration April 2022

SIZE	1 gallon (3.8L)
DECON #	8811

Characteristic	Specifications	Test Result
Appearance	Clear colorless to light yellow liquid	Pass
Odor	Nectar	Pass
ortho-Phenylphenol:	3.23 - 3.57%	3.33
ortho-Benzyl-para-chlorophenol	2.88 - 3.18%	2.94
pH: Concentrate	2.00 - 2.40	2.15
pH: 1:128 Dilution	2.50 - 2.90	2.81
Density (lb/gal):	8.74 - 8.82	8.76

Reviewed and Approved by:

Date: April 2020

A handwritten signature in blue ink that reads 'Karen Wend'.

Karen Wend
Manager of Quality and Technical Services

BSE/TSE Statement: There are no materials of animal origin in the raw materials or in the process aides used to manufacture any Decon products.