




K:\GCA\223686

THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION AND IS A COPYRIGHT OF WHEATON INDUSTRIES INC. IT MAY NOT BE DUPLICATED FOR ANY PURPOSE IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN CONSENT OF WHEATON INDUSTRIES INC. ALL RIGHTS RESERVED.©
ELECTRONIC OR PAPER COPIES OF THIS PRINT OUTSIDE OF WHEATON'S SYSTEM ARE UNCONTROLLED AND ARE FOR REFERENCE ONLY UNLESS OTHERWISE STAMPED.

THIRD ANGLE PROJECTION



ISOMETRIC VIEW



MATERIAL: TYPE I CLEAR BOROSILICATE GLASS

3. VIAL IS SUITABLE FOR FREEZE DRYING.
2. INCH DIMENSIONS ARE FOR REFERENCE ONLY.
1. DIMENSIONS ARE IN MILLIMETERS & [INCHES].

----NOTES----

TITLE	10 mL VIAL, SERUM, TYPE I CLR GLASS
CAPACITY	10cc (MEAN 15.1cc OVERFLOW)
FINISH	G.P.I. 2710-20A W/ BLOWBACK
SCALE: FULL UNLESS NOTED. DO NOT SCALE THIS DRAWING.	

DATE DRAWN	02/25/03	PROJ#	
DRAWN BY	B. LLOYD	APPROX WGT	
CHECKED BY	TRJ 02/28/03	AUTHORIZED	SS 02/28/03
WHEATON® 1501 North 10th Street, Millville, NJ 08332-2038			

-	REDRAWN IN SOLIDWORKS UPDATED TO W.I. STANDARDS	08/16/16 GT
B	2013011411102976	01/17/13 TP
A	395516747972	02/25/03 BL
REV	ECN #	DATE
SIZE	DWG.NO: 223686	
B		