

Product Description	:	Beaker - Polypropylene, Low Form, Screen Printed
Capacity	:	10,000 ml
Product Code	:	CH0139H
Lot Number	:	BA 090 34
Country of Origin	:	Manufactured in India

Compliance: This is to certify that the product is manufactured in ISO 9001 : 2015 certified plant. Products are inspected and controlled in accordance with current applicable product specifications. Inspection records are reviewed and maintained.

Material	:	Polypropylene - Comply with the FDA: CFR Title 21, 177.1520 olefin polymer, product complies with the requirements of Regulation EU/10/2011 for plastics intended to come to contact with foodstuff and the requirement of Overall Migration Limit (OML) of 60 mg/kg as per mentioned in EU/10/2011 and IS 9845: The material used conforms to amendments to 94/62/EC (Amendment number EU/2015/720, 2004/12/EC, 2005/20/EC, 219/2009). This grade meets the requirement of less than 100 ppm for total incidental Cadmium, Lead, Chromium, and Mercury, we do not intentionally use or add phthalates, plasticizers, fillers, or flexibilisers.
BSE/TSE	:	The above product is manufactured without using any raw material produced or derived from animal origin.
DNase & RNase Free	:	Not Applicable
Cytotoxicity	:	Not Applicable
Latex	:	No Latex was used in the manufacturing of the product. This also includes all the packaging and shipping materials used in the production of the items.
BPA Statement	:	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
Storage Conditions	:	Store in cool and dry location away from direct sunlight and heat in normal warehouse conditions.

(An electronic generated document does not require any signatures)