



Certificate of Conformity

VWR WATCH GLASSES, PTFE

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat. No.	Ø ext. (mm)	Pk
216-0454	20	4
216-0155	30	1
216-0156	40	1
216-0157	50	1
216-0158	65	1
216-0159	75	1
216-0160	80	1
216-0161	100	1
216-0162	125	1
216-0163	150	1
216-0455	200	1

Material: PTFE conforming to FDA and USP Class VI requirements

Properties: Good chemical resistance

Temperature Resistance: -200° up to +280°C

Compliant to: FDA

PTFE used in the manufacture of these products are pure, virgin material which conforms to FDA and USP Class VI requirements. It does not contain, or use in its manufacture any Class 1 ODC according to the Clean Air Act 1990, PFOA, POP, DEHIP and BPA, no animal material is used in the manufacture or is present in this material. PTFE used conforms to current REACH requirements.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **October 2017**