



Product Specification Sheet

VWR ADHESIVE FILM CLEANING MATS FOR CLEANROOMS, 30-LAYER, PURESTEP

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	L×W mm	Thickness mm	Colour	Pk
115-0061	457×1143	1,8	Blue	8
115-0062	457×1143	1,8	Grey	8
115-0063	457×1143	1,8	White	8
115-0055	457×914	1,8	Blue	8
115-0056	457×914	1,8	Grey	8
115-0057	457×914	1,8	White	8
115-0047	610×762	1,8	Blue	4
115-0048	610×762	1,8	Clear	4
115-0049	610×762	1,8	Grey	4
115-0050	610×762	1,8	White	4
115-0067	610×914	1,8	Blue	8
115-0068	610×914	1,8	Grey	8
115-0069	610×914	1,8	White	8
115-0073	660×1143	1,8	Blue	8
115-0074	660×1143	1,8	Grey	8
115-0075	660×1143	1,8	White	8
115-0085	914×1143	1,8	Blue	8



Product Specification Sheet

VWR ADHESIVE FILM CLEANING MATS FOR CLEANROOMS, 30-LAYER, PURESTEP

Cat. No.	L×W mm	Thickness mm	Colour	Pk
115-0086	914×1143	1,8	Grey	8
115-0087	914×1143	1,8	White	8
115-0091	914×1524	1,8	Blue	4
115-0092	914×1524	1,8	Grey	4
115-0093	914×1524	1,8	White	4
115-0079	914×914	1,8	Blue	8
115-0080	914×914	1,8	Grey	8
115-0081	914×914	1,8	White	8

Other Specifications:

Film : Low density Polyethylene
File Thickness/Layer : 1.8 mil / 50 micron
Adhesive : Acrylic
Adhesive Thickness : 0.3 mil/8 micron
Peel Adhesion : 6.5 oz./180 g/25mm $\pm$ 20 (individual layers)
24 oz./680 g/25mm $\pm$ 50 (bottom layers)
Tensily Strength : >1.2 kg/cm
Heat Resistance : 158° F/70° C -48 hours

Other Characteristics:

- Anti-microbial, destroys germs and fungi
- Non-slip, colour-coded
- Reduced maintenance and cleaning costs

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt, ask us.