

Product Specification

Test	Specification
Assay (by GC)	>= 95 %
Residue after Evaporation	<= 2 ppm
Substances Darkened by H2SO4	passes test
Water (H2O)	<= 0.01 %
Ultraviolet Absorbance,(1.00-cm path vs water; curve smooth throughout stated range with no extraneous impurity peaks).	
at 210 nm	<= 0.2
at 220 nm	<= 0.08
at 254 nm	<= 0.008
at 280 nm	<= 0.005
at 350 nm	<= 0.005
UV Cut-off, nm	<= 192

For Laboratory, Research or Manufacturing use.
For use in High Performance Liquid Chromatography.
Filtered through a 0.2 micron filter.Packaged under Nitrogen.
Because of the possible presence of aluminium particles it is advised to filter before use, over a 0.2 or 0.45 µm filter.



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000