


Recombinant Human UCHL1 Protein (A270577)

Specifications:

Name:	Recombinant Human UCHL1 Protein
Applications:	ELISA, MSD, Luminex, Simoa assays
Expression System:	Escherichia coli
Nature:	Recombinant
Protein Species:	Human
Protein Length:	Full length protein.
Sequence:	The vector used in the recombinant expression of this protein adds about 5 kDa to the molecule, including an N-terminal His-tag.
Tag:	His Tag (N-terminus)
Product Form:	Liquid
Formulation:	Supplied in 10 mM Phosphate Buffer, pH 7.5, with 6 M Urea.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

Images:



Coomassie brilliant blue stained SDS-PAGE gel of Recombinant Human UCHL1 Protein (A270577). The lanes contain: [Lane 1] Recombinant Human UCHL1 Protein (A270577) and [Lane 2] protein standards.