

Certificate of Conformity

VWR Catalogue Number	211-0255
Description	Cryoboxes, for light-sensitive specimens, 100 compartments, for Microcentrifuge tubes 0,5/1,0 and 2,0 ml, 140 mm width, 140 mm depth, 60 mm height, Black
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2022-11-04

Quality System Compliance

Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polypropylene.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.

Latex Statement Not Applicable.

BPA Statement Not Applicable.

DEHP Statement Not Applicable.

RoHS No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.

REACH Statement No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.

Storage Conditions Indefinite.

Shelf Life Indefinite.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.