



Certificate of Conformity

VWR Mega Microtube Rack

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat. No.	For Tubes	No of holes	WxDxH	Colour	Pk
211-0364	0,2/0,5/1,5/2,0/5,0ml	216	337×267×91 mm	White	1
211-0363	0,2/0,5/1,5/2,0/5,0ml	216	337×267×91 mm	Blue	1
212-0661	0,2/0,5/1,5/2,0/5,0ml	432	673×230×80 mm	White	1
212-0660	0,2/0,5/1,5/2,0/5,0ml	432	673×230×80 mm	Blue	1

Material: Polypropylene

Other Characteristics:

- For tubes of diameter 10 to 13 mm
- High quantity sample storage eliminates the need for multiple racks
- Alphanumeric grid for quick location and retrieval of tubes
- Sturdy, lightweight material for easy transport

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **February 2019**