

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Low Form Griffin Beakers
Country of Origin	Manufactured in Czech Republic
Date of Issue (yyyy-mm-dd)	2025-07-28

VWR® Low form griffin beakers

VWR Catalogue Number	Capacity	Height	Ext.Ø	Graduation Range	Unit
10754-696	10 ml (0.34 oz.)	34 mm (1.34")	25 mm	—	Pack of 12
					Case of 48
10754-698	20 ml (0.68 oz.)	41 mm (1.6")	32 mm	5 - 15 ml	Pack of 12
					Case of 48
10754-700	30 ml (1 oz.)	53 mm (2.1")	35 mm	5 - 25 ml	Pack of 12
					Case of 48
10754-946	50 ml (1.7 oz.)	53 mm (2.1")	42 mm	10 - 40 ml	Pack of 12
					Case of 48
10754-948	100 ml (3.4 oz.)	70 mm (2.76")	50 mm	20 - 80 ml	Pack of 12
					Case of 48
10754-950	150 ml (5 oz.)	81 mm (3.19")	57.5 mm	20 - 140 ml	Pack of 12
					Case of 48
10754-952	250 ml (8.5 oz.)	88 mm (3.46")	68 mm	25 - 200 ml	Pack of 12
					Case of 48
10754-954	400 ml (14 oz.)	106 mm (4.17")	77 mm	25 - 325 ml	Pack of 12
					Case of 48
10754-956	600 ml (20 oz.)	122 mm (4.8")	88 mm	50 - 500 ml	Pack of 6
					Case of 36
10754-958	800 ml (27 oz.)	136 mm (5.35")	99 mm	50 - 750 ml	Pack of 6
					Case of 24
10754-960	1000 ml (34 oz.)	145 mm (5.7")	107 mm	100 - 1000 ml	Pack of 6
					Case of 24
10754-962	1500 ml (50.7 oz.)	165 mm (6½")	120 mm	200 - 1400 ml	Pack of 4
					Case of 16
10754-760	2000 ml (68 oz.)	188 mm (7.4")	131 mm	200 - 2000 ml	Pack of 4
					Case of 8
10754-762	4000 ml (135.3 oz.)	250 mm (9.8")	160 mm	500 - 3500 ml	Pack of 1
					Case of 6

VWR® Low form griffin beakers sample kit

VWR Catalogue Number	Capacity	Unit
75870-432	25 ml, 50 ml, 125 ml, 250 ml, 300 ml	Kit of 5

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001 and ISO 28000** standard.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Test 1: Hydrolytic resistance of glass grains at 98°C

Test reference: ISO 719:2020
 Acceptance Criteria: HGB1: < 0.10 ml/g glass
 Test Result: 0.02 ml/g glass
 Compliance: YES

Test 2: Volume and Dimensional verification

Test reference: Technical drawings
 Acceptance Criteria: Technical drawings
 Test Result: OK
 Compliance: YES

The product is in conformance with the requirements of the following standards and regulations:

- **Glass characteristics:** ISO 3585 Borosilicate glass 3.3 –Properties.
- **Technical standards for products:** ISO 3819 Laboratory glassware - Beakers
 ISO 4794 Laboratory glassware - Methods for assessing the chemical resistance of enamels used for colour coding and colour marking
- **Glass containers for pharmaceutical use:** European Pharmacopoeia 10th, Edition Glass Type I & US Pharmacopoeia-Glass Type I.

Product Specifications

Material	SIMAX® Borosilicate glass. The FDA classifies glass in general as GRAS (Generally Recognized as Safe) according to Federal Register/Vol 48, No. 27/Feb. 08.1983/page 5718, Rules and Regulations.
Sterilization	The products are non-sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Store at room temperature in normal warehouse conditions; Protect from any possible contamination.
Shelf-Life	Indefinite.

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	July 2025