

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Disposable Transfer Pipettes
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2026-03-31

VWR Catalogue Number	Description	Capacity [mL]	Overall length	Sterility	Drops per [mL]	Bulb draw [mL]	Pk
16001-170	Short thin stem	1.2	6.4 cm (2 1/2")	Non sterile	25	0.9	500/box
16001-192	Extended fine tip, small bulb	1.5	10.4 cm (4")	Non sterile	50	1.0	400/box
414004-018	Extended fine tip, small bulb	1.5	10.4 cm (4")	Sterile	50	1.0	1/bag, 400/box
414004-019	Extended fine tip, small bulb	1.5	10.4 cm (4")	Sterile	50	1.0	20/bag, 500/box
16001-174	General-purpose	1.7	8.7 cm (3 7/16")	Non sterile	25	0.9	500/box
414004-030	General-purpose	1.7	8.7 cm (3 7/16")	Sterile	25	0.9	1/bag, 500/box
414004-031	General-purpose	1.7	8.7 cm (3 7/16")	Sterile	25	0.9	20/bag, 500/box
414004-007	Narrow stem, large bulb	4.0	8.5 cm (3 3/8")	Non sterile	25	3.0	500/box
414004-033	Blood bank	4.0	12.9 cm (5")	Non sterile	25	1.9	400/box
414004-038	Blood bank	4.0	12.9 cm (5")	Sterile	25	1.9	20/bag, 500/box
414004-039	Blood bank	4.0	12.9 cm (5")	Sterile	25	1.9	1/bag, 400/box
16001-178	Standard, thin stem	4.6	15.5 cm (6 1/8")	Non sterile	22	3.3	500/box
414004-008	Narrow stem, standard	4.6	15.5 cm (6 1/8")	Sterile	22	3.3	1/bag, 500/box
414004-009	Narrow stem, standard	4.6	15.5 cm (6 1/8")	Sterile	22	3.3	5/bag, 500/box
414004-010	Narrow stem, standard	4.6	15.5 cm (6 1/8")	Sterile	22	3.3	10/bag, 500/box
414004-011	Narrow stem, standard	4.6	15.5 cm (6 1/8")	Sterile	22	3.3	20/bag, 500/box

16001-194	Fine tip, standard	5.0	14.9 cm (5 ⁷ / ₈ "	Non sterile	50	3.3	500/box
414004-001	Blood bank	5.0	15.5 cm (6 ¹ / ₈ "	Non sterile	25	1.9	500/box
414004-002	Blood bank	5.0	15.5 cm (6 ¹ / ₈ "	Sterile	25	1.9	1/bag, 500/box
414004-003	Blood bank	5.0	15.5 cm (6 ¹ / ₈ "	Sterile	25	1.9	20/bag, 500/box
414004-023	Extended fine tip, large bulb	5.0	15.2 cm (6"	Non sterile	50	3.3	250/box
414004-042	Extended fine tip, large bulb	5.0	15.2 cm (6"	Sterile	50	3.3	20/bag, 500/box
16001-186	Long	5.5	23 cm (9"	Non sterile	25	2.3	500/box
414004-026	Long	5.5	23 cm (9"	Sterile	25	2.3	1/bag, 400/box
414004-028	Long	5.5	23 cm (9"	Sterile	25	2.3	20/bag, 400/box
414004-020	Fine tip, standard	5.8	14.7 cm (5 ²⁵ / ₃₂ "	Non sterile	50	3.3	500/box
414004-021	Fine tip, standard	5.8	14.7 cm (5 ²⁵ / ₃₂ "	Sterile	50	3.3	1/bag, 500/box
414004-022	Fine tip, standard	5.8	14.7 cm (5 ²⁵ / ₃₂ "	Sterile	50	3.3	20/bag, 500/box
414004-048	Fine tip, standard	5.8	14.7 cm (5 ²⁵ / ₃₂ "	Non sterile	50	3.3	10/bag, 500/box
414004-004	Standard	7.5	15.5 cm (6 ¹ / ₈ "	Non sterile	20	3.2	500/box
414004 005	Standard	7.5	15.5 cm (6 ¹ / ₈ "	Sterile	20	3.2	1/bag, 500/box
414004-006	Standard	7.5	15.5 cm (6 ¹ / ₈ "	Sterile	20	3.2	20/bag, 500/box
16001-190	General-purpose, large bulb	8.0	15.2 cm (6"	Non sterile	20	4.6	400/box
414004-046	General-purpose, large bulb	8.0	15.2 cm (6"	Sterile	20	4.6	20/bag, 400/box
414004-047	General-purpose, large bulb	8.0	15.2 cm (6"	Sterile	20	4.6	1/bag, 400/box
414004-032	Large aperture	9.3	13.7 cm (5 ⁷ / ₁₆ "	Non sterile	18	4.8	250/box
414004-012	Narrow stem, large bulb	15.0	15.5 cm (6 ¹ / ₈ "	Non sterile	19	8.0	250/box
414004-041	General-purpose, large bulb	15.0	15.5 cm (6 ¹ / ₈ "	Sterile	19	8.0	1/bag, 250/box
16001-196	Extra long	23.0	30 cm (11 ¹³ / ₁₆ "	Non sterile	22	7.3	100/box

414004-027	Extra long	23.0	30 cm (11 ¹³ / ₁₆ "	Sterile	20	7.3	1/bag, 100/box
414004-029	Extra long	23.0	30 cm (11 ¹³ / ₁₆ "	Sterile	20	7.3	20/bag, 100/box

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with the [ISO 13485:2016](#) standard.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Low Density PolyEthylene (LDPE) DJ200A (Virgin).
Sterilization	The product labelled as sterile is EB (Electron Beam) irradiated and meets a minimum requirement of 10 ⁻⁶ SAL with a specified dose range of 15 - 40 kGy.(which is not applicable for non-sterile products).
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin, and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.

REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Avoid high temperature & moisture. No heavily loading during transport or storage.
Shelf-Life	4 years. .(which is not applicable for non-sterile products).

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	October 2021
2	2.0	Modified Quality system compliance, disclaimer, REACH, shelf life, removed FDA	April 2026